

HBT ELISA TEST KIT FOR HUMAN IFN-gamma 2 x 96 DETERMINATIONS

The Hbt Human IFN-gamma kits have been developed for the quantitative measurement of natural and recombinant Human IFN-gamma in serum, plasma and culture medium.

PRINCIPLE OF THE TEST

The Hbt Human IFN-gamma ELISA test kit is a solid-phase enzyme-linked immunosorbent assay based on the sandwich principle. Samples and standards are incubated in microtiter wells coated with antibodies recognizing Human IFN-gamma. During this incubation Human IFN-gamma is captured by solid bound antibody. Unbound material present in the sample is removed by washing. Next biotinylated second antibody (tracer) to Human IFN-gamma is added to the wells. If IFN-gamma was present in the sample, the tracer antibodies will bind to the captured IFN-gamma. Excess tracer is removed by washing. Next a streptavidin-peroxidase conjugate is applied to the wells, this conjugate reacts specifically with the biotinylated tracer antibody bound onto the captured IFN-gamma. Excess streptavidin-peroxidase conjugate is removed by washing and the substrate, tetramethylbenzidine (TMB) is added to the wells. Colour develops proportionally to the amount of human IFN-gamma present in the sample. The enzyme reaction is stopped by the addition of citric acid and the absorption at 450 nm is measured with a spectrophotometer. A standard curve is obtained by plotting the absorptions versus the corresponding concentrations of the known standards. The human IFN-gamma concentration of samples with unknown concentrations, which are run concurrently with the standards, can be determined from the standard curve.

SPECIAL FEATURES OF THE KIT

- Ready-to-use (i.e. pre-coated microwells).
- High specificity.
- High reproducibility
- High sensitivity. The minimum concentration which can be measured is 15 pg/ml of IFN-gamma.
- Large measurable concentration range. Standard curve from 20-5,000 pg/ml.
- Efficient format. Two plates with each twelve 8-well strips allow free choice of batch size for the assay.
- Simple, rapid procedure. Four pipetting steps are required to complete the assay. Working time 4½ hours.

AVAILABILITY

The Hbt Human IFN-gamma tests are available in kits for 2 x 96 and 10 x 96 determinations.

PRODUCT NUMBER: HK309 Hbt Human IFN-gamma ELISA Kit, 2 x 96 determinations

Also available:

PRODUCT NUMBER: HK030m Hbt Human IFN-gamma ELISA Kit (Medium), 1 x 96 determinations

PRODUCT NUMBER: HK030s Hbt Human IFN-gamma ELISA Kit (Serum), 1 x 96 determinations

PRODUCT NUMBER: DK1309 Hbt Human IFN-gamma DECA Kit, complete ELISA reagent set for 10 x 96 determinations.

For research purposes only.

Caution: Not for use in humans.